



CAH (21-HYDROXYLASE DEF)

14755X

CPT Code(s): PROF

Preferred Specimen(s)

****PLEASE NOTE THIS TEST IS FOR NON NY CLIENTS. NY CLIENTS PLEASE ORDER SOC 15063X****

5 mL whole blood from TWO EDTA lavender top tubes.

Minimum Volume

3 mL whole blood from a EDTA lavender top tube

Alternate Specimen(s)

5 mL whole blood from any of the following: ACD solution B yellow top tube, EDTA royal blue top tube, or sodium heparin green top , or ACD solution A yellow top tube.

For tissue biopsy samples please contact client services.

Please contact client services prior to submitting 100 ng Extracted DNA.

Transport Container

5 mL whole blood from 2 EDTA lavender top tubes

Transport Temperature

Room temperature

Specimen Stability

Room temp: 8 days

Refrigerated: 8 days

Frozen: Unacceptable

Methodology

PCR, electrophoresis, and minisequencing

Setup Days

Tuesday

Report Available

8 days

(The CPT codes provided are based on AMA guidelines and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.)